- DE** Bedienungsanleitung
- GB** Instruction Manual
- FR** Mode d'emploi
- ES** Instrucciones de uso

B BRESSER®

XPD-101

Expeditionsmikroskop

Art. No. 5790500

VORSICHT!

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All parts:

- 1 Eyepiece
- 2 Eyepiece socket
- 3 objective revolver
- 4 On/Off switch
- 5 Object stage
- 6 Object stage-Adjustment
- 7 Focussing wheel
- 8 Battery compartment cover
- 9 Fixing screws for battery compartment cover
- 10 Object stage with object plate
- 11 Object fixing clamp

1. Choosing a place for the microscope

Before starting, choose a sufficient location for microscoping. Please take care that your microscope stands on an even and solid ground.

2. Getting started

The expedition microscope is preassembled and ready to use. Before use, just insert the supplied 10x eyepiece in the eyepiece holder and 2 batteries (type AA, 1.5V) in the battery compartment. Open the battery compartment cover by loosening the two knurled screws and slide the cover plate away to the side. Insert the batteries as required and close the cover again.

Pull the eyepiece holder out and lock it in the upper position by a $\frac{1}{4}$ turn counterclockwise. This is the only way to reach proper focus, otherwise the travel of the focus mechanics will not be sufficient. Insert the eyepiece into the eyepiece socket. Insert a microscope slide into the clamp of the stage. The specimen with cover glass is inserted upright (i.e. cover glass on top). The objectives are suitable for slides with thickness between 1.0-1.5mm. No special plankton chambers etc. with cover glass bottom are required, thus reducing material cost significantly.

3. Observation

Switch on the light with the toggle switch. Start with the lowest magnification (4x objective) to scan through your sample. The objectives 4x and 10x produce an image with a type of dark field contrast, i.e. sample components appear bright against a dark background. The 40x objective is equipped with a phase ring so that samples are viewed in phase contrast. Objects appear brightly contrasted against a grey background. In phase contrast, even the smallest, colourless and translucent objects and organelles (e.g. flagella) can be detected more easily than in brightfield. For plankton this is the method of choice.

To store or transport the microscope in the smallest space possible, the microscope can be disassembled without tools. Loosen the knurled screws that connect the illumination housing with the base, and then loosen the knurled screws on the cross table and remove the crosstable. Store microscope body and parts in a padded bag or case.

If the illumination in the field is not ready (e.g. batteries empty), the microscope can also be operated with ambient light only. To do this, simply remove the housing of the illumination. The

daylight entering from above is sufficient for working with low magnifications. When used like this, the microscope gives a normal brightfield image.

The microscope can be used hand-held when observing for a short time. If the microscope is to be used more often or for a longer period of time, the device can be mounted on a photo tripod or small table tripod for convenient viewing, or you can mount it on any solid surface with a table clamp. For this purpose, there is a ¼" camera tripod thread in the base plate.

4. Cleaning and maintenance

Your microscope is a high-quality optical instrument. Please take care that neither dust nor fingerprints take place on the optical surfaces. If necessary, clean the microscope body with a slightly moistened soft cloth. Use a separate clean cloth with some drops of alcohol to clean the optical surfaces. After using, the microscope should be protected with its dust cover.

Note:

Please consider that a sufficient treated microscope keeps its value and optical quality for many years.

5. EC Declaration of Conformity

 Bresser GmbH has issued a „Declaration of Conformity“ in accordance with applicable guidelines and corresponding standards. The full text of the EU declaration of conformity is available at the following internet address: www.bresser.de/download/5790500/CE/5790500_CE.pdf

6. Only for EU countries

Do not dispose of electric equipment together with household waste material!

 In observance of European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE) and its implementation in accordance with national law, electric equipment that have reached the end of its life must be collected separately and returned to an environmentally compatible recycling facility.

Discharged batteries and damaged re-chargeable batteries must be disposed of at special battery collection points. Information is available from your local disposal agent or local authority regarding the disposal of devices or batteries.

7. Warranty & Service

The regular guarantee period is 10 years and begins on the day of purchase. You can consult the full guarantee terms and details of our services at www.bresser.de/warranty_terms.

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