

Art. No. 88-54500

- DE** Bedienungsanleitung
- GB** Operating Instructions
- FR** Mode d'emploi
- NL** Handleiding
- IT** Istruzioni per l'uso
- ES** Instrucciones de uso

DE	Bedienungsanleitung	4
GB	Operating Instructions	11
FR	Mode d'emploi	18
NL	Handleiding	25
IT	Istruzioni per l'uso	32
ES	Instrucciones de uso	39

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DANGER to your child!

This device contains electronic components that are powered by either a mains connection or batteries. Never leave a child unsupervised with this device. The device should only be used as per these instructions otherwise there is a serious RISK of ELECTRICAL SHOCK.

Children should only use this device under supervision. Keep packaging materials (plastic bags, rubber bands, etc.) away from children. There is a risk of SUFFOCATION.

FIRE-/ DANGER OF EXPLOSION!

Do not expose the device to high temperatures. Use only the mains adapter supplied or those battery types recommended. Never short circuit the device or batteries or throw into a fire. Exposure to high temperatures or misuse of the device can lead to short circuits, fire or even explosion!

Do not subject the device to temperatures exceeding 60 C.

RISK of material damage.

Never take the device apart. Please consult your dealer if there are any defects. The dealer will contact our service centre and send the device in for repair if needed.

Cleaning instructions

Remove the device from its energy source before cleaning (remove plug from socket / remove batteries)

Clean the exterior of device with a dry cloth. Do not use cleaning fluids so as to avoid causing damage to electronic components.

Protect the device from dust and moisture. Store the device in the bag supplied or in its original packaging. Batteries should be removed from the device if it is not going to be used for a long period of time.

DISPOSAL

Dispose of the packaging material/s as legally required. Consult the local authority on the matter if necessary.

Do not dispose of electrical equipment in your ordinary refuse. The European guideline 2002/96/EU on Electronic and Electrical Equipment Waste and relevant laws applying to it require such used equipment to be separately collected and recycled in an environment-friendly manner.

Empty batteries and accumulators must be disposed of separately. Information on disposing of all such equipment made after 01 June 2006 can be obtained from your local authority.

Declaration of Conformity

CE Meade Instruments Europe GmbH & Co KG, resident in 46414 Rhede/Westf., Gutenbergstr. 2, Germany, explains the agreement with in the following specified guidelines for this product:

EN 55022:2006+A1:2007

EN 61000-3-2:2006

EN 61000-3-3:2008

EN 55024:1998+A2:2003

Product description:
USB Universal Microscope
Model /Description:
Digital Microscope

Rhede, 08-01-2010

Meade Instruments Europe
GmbH & Co. KG

Helmut Ebbert
Managing director

All parts:

- 1 Camera shutter
- 2 Reduce brightness (-)
- 3 Increase brightness (+)
- 4 Set light mode / light off
- 5 Focus Ring
- 6 Attachment Piece
- 7 Illuminator (12 LED Lights)
- 8 Lens
- 9 USB Cable
- 10 USB Connector
- 11 Stand
- 12 Focusing knob
- 13 Illuminated table
- 14 Preparation holders
- 15 Box with Object holder,
glass covers and Adhesive labels

General

This is a digital reflected light microscope. You hold it in your hand and can place the bottom section (attachment piece) on all kinds of things in order to look at them. Observe, for example, leaves, microorganisms, your skin or hair, and much more. It

works best when the thing that you're looking at (also called the "object") is flat. You can also view the enlarged pictures on your computer, as well as take and save them there.

Installation

Insert the product CD into your PC's DVD/CD drive. The driver installation starts automatically. Plug the hand microscope's USB connector (10) into your computer's USB slot. The lighting (7) turns on and your PC detects the hardware, which is then installed. Soon, the "AMCAP" icon appears on the desktop. Now you can use the hand microscope.

Live Observation

Press the camera shutter (1) for your hand microscope. A (in general blurry) live image is displayed on the monitor.

Hold the hand microscope by the casing and place the attachment

piece (6) on an object, for example a piece of paper with writing on it. Turn the focus ring (5) to make the live picture sharper (this is called focusing). For a flat object, there are two focus settings with sharp images, which correspond to two different magnifications. For low-power magnification, the lens (8) is positioned high, away from the object. For high-power magnification, it is positioned lower, closer to the object. You can adjust the magnification from low to high by turning the focus ring clockwise. To turn it from low magnification to high, turn the focus ring counter-clockwise. You'll only know when you have the exact measurement value when you've achieved a clear picture of an object (e.g., as shown on your computer screen or printed out on a piece of paper).

Turn the microscope until you have a picture that is straight and right side up.

Light Mode Settings:

You can select from 4 different light settings with “MODE”:

- white light
- white and red light
- white and yellow light
- white, red and yellow light

The light can be shut off altogether using the “MODE” key as well. Regulate the brightness using the “+” and “-” keys on the device. (Press and hold the keys!)

Taking Pictures

Using the camera shutter (1), you can take a picture and save it as a BMP file.

1. Press the camera shutter (1).
2. The “SnapShotView” window appears on the screen with a picture.
3. To save the image, click “File” and “Save”.

Making Movies

The “AMCAP” program allows you to make movies with the hand microscope and save them as AVI files.

1. Click on “File” and “Set Capture File...”; specify the name of the AVI file with the file extension “.avi”. For example: “experiment1.avi”.
2. The “Set File Size” window appears on the screen. Here, specify the maximum file size.
3. You can prepare to film with “Capture” and “Start Capture” in the menu.
4. Start filming with “OK” in the “Ready to Capture” window.
5. Under the “Capture” menu, you can end your filming with “Stop Capture”.
6. If you want to record a new film, follow Step 1 and specify a new AVI file with a new name. Otherwise, the file will overwrite your film.
7. You can watch your film using a playback program for multimedia files.

Magnifications

In the lower magnification, a picture includes about 10,5 mm x 14 mm of the object. The higher magnification includes about 1 mm x 1,4 mm. In this way, the higher magnification is about five times stronger than the lower one. When, for example, you print a picture that is 28 cm wide on a piece of paper, the magnification is about 20x (low) or 200x (high).

Deactivation and Storage

Close the “AMCAP” window on your PC screen. Now, you can remove the USB connector (10) from your computer’s USB port. You can store your hand microscope in the storage case until the next time you want to use it. This will protect it from dust.

Technical Information

- Digital hand microscope with computer connection (USB)
- Magnification: 20x & 200x
- Bright illumination via 12 LEDs
- Power supply via USB
- Image preview:
 - 15 fps (USGA: 1.280x1.024)
 - 30 fps (GUXGA: 800x600)
- Size: 54x54x104 mm
- Weight: 144g

System Requirements

Windows XP with Service Pack 3 (on CD-ROM), Windows Vista, Windows 7 - with DirectX 9.x (on CD-ROM), a minimum of 1 GB RAM, free USB 2.0 port.

Photomizer SE Software

Photomizer SE Software can be downloaded free of charge from:
<http://www.bresser.de/downloads/support/software/photomizer.zip>

Experiments with the Handheld Digital Microscope

Experiment No. 1: Black and White Print

Objects:

1. a small piece of paper from a newspaper with a black and white picture and some text,
2. a similar piece of paper from a magazine.

Place both pieces of paper next to each other on a table. Set your microscope to the lowest magnification and place it on the pieces of paper, first on the newspaper and then on the magazine.

Compare: The letters on the newspaper look frayed and broken, since they are printed on raw, low-quality paper. The letters on the magazine look smoother and more complete. The pictures in the newspaper are made up of many tiny dots, which appear slightly smudgy. The pixels (half-

tone dots) of the magazine picture are clearly defined.

Experiment No. 2: Color Print

Objects:

1. a small piece of color printed newspaper,
2. a similar piece of paper from a magazine.

Place both pieces of paper next to each other on a table. Set your microscope to the lowest magnification and place it on the pieces of paper, first on the newspaper and then on the magazine.

Compare: The colored pixels of the newspaper often overlap. Sometimes, you'll even notice two colors in one pixel. In the magazine, the dots appear clear and rich in contrast. Look at the different sizes of the pixels.

Experiment No. 3: Textile fibers

Objects and accessories:

1. threads from various fabrics (e.g. cotton, linen, sheep's wool, silk, rayon, etc.),
2. two needles.

Place the different threads on a table and use the needles to fray them a bit. Dampen the threads with a little water. Set your microscope to the lowest magnification and place it on the threads, one at a time.

Compare: Cotton fibers come from a plant, and look like a flat, twisted ribbon under the microscope. The fibers are thicker and rounder at the edges than in the middle. Cotton fibers are basically long, collapsed tubes. Linen fibers also come from a plant, and they are round and run in one direction. The fibers shine like silk and exhibit countless bulges on the thread. Silk comes from an animal and is made up of solid fibers that are small in diameter, in contrast to

the hollow plant-based fibers. Each fiber is smooth and even and looks like a tiny glass tube. The fibers of the sheep's wool also come from an animal. The surface is made of overlapping sleeves that look broken and wavy. If possible, compare sheep's wool from different weaving mills. In doing so, take a look at the different appearance of the fibers. Experts can determine which country the wool came from by doing this. Rayon is a synthetic material that is produced by a long chemical process. All the fibers have solid, dark lines on the smooth, shiny surface. After drying, the fibers curl into the same position. Observe the differences and the similarities.

Experiment No. 4: Table Salt

Object: normal table salt.

Place a sheet of black paper on a desk. Sprinkle a few grains of salt on the paper and place the microscope

on top of them. Look at the salt crystals using the lowest magnification of your microscope.

Observe: The crystals look like tiny dice and all have the same shape.

Experiment No. 5: Leaves and Needles

Object: 3-4 different leaves or needles from deciduous trees or fir trees.

When you go for a walk in the forest with your parents, you can collect different types of leaves and needles.

At home, place them next to each other on a white sheet of paper. Place your microscope on top of them and look at the different leaves and needles with the lowest magnification.

Observe: The leaves of the deciduous trees have different but more or less regular sections that are separated by lines. These are called "cells." Most often, the underside of the leaf looks different than the top, and the

color is brighter. The stalk of the leaf runs through the middle. At its thicker end, there is a “lump” with a bulge. That is the part that connected the leaf to the tree, before it fell away. Some leaves also have a stalk upon which multiple leaves grow from other stalks.

Fir needles are long, thin and round. Like the leaves of deciduous trees, they have a light bulge on one side, where they grew from the tree. They do not have individual “cells,” however, but look like they grew in one part. However, when you look more closely, you can see that the needle has many sections. These sections come from the step-by-step growth of the needles.

In this way, you can look at many more objects, such as small organisms (flies, spiders, etc.) or other things from your daily life. Simply put everything on a flat surface (a desk) and place the microscope on top.

Or have you already looked at the hair on your head? No? Then run the hand microscope through your hair. It's quite funny and surprising, what can be hidden in there.

You can discover so many things that you did not know before. Just give it a try!

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