Biolux

Mikroskop

Art. No. 9810302

- DE Bedienungsanleitung
- EN Operating instructions

Fig. 1










Fig. 2

SERVICE AND WARRANTY:


www.bresser.de/warranty_terms

MANUAL DOWNLOAD:

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
MICROSCOPE GUIDE:

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EXPERIMENTS:

www.bresser.de/downloads

MICROSCOPE FAQ:

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DE	Bedienungsanleitung.....	5
EN	Operating instructions	10

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Operating Instructions

WARNING! Not for children under three years. Choking hazard – small parts. Contains functional sharp edges and points!

ATTENTION: Only suitable for children of at least 8 years of age. Instructions for parents or other responsible persons are enclosed and must be followed. Keep packaging as it contains important information.

General Warnings

- **Choking hazard** – This product contains small parts that could be swallowed by children. This poses a choking hazard.
 - **Risk of electric shock** – This device contains electronic components that operate via a power source (batteries). Only use the device as described in the manual, otherwise you run the risk of an electric shock.
 - **Risk of fire/explosion** – Do not expose the device to high temperatures. Use only the recommended batteries. Do not short-circuit the device or batteries, or throw them into a fire. Excessive heat or improper handling could trigger a short-circuit, a fire or an explosion.
 - **Risk of chemical burn** – Make sure you insert the batteries correctly. Empty or damaged batteries could cause burns if they come into contact with the skin. If necessary, wear adequate gloves for protection.
- Leaking battery acid can lead to chemical burns. Avoid contact of battery acid with skin, eyes and mucous membranes. In the event of contact, rinse the affected region immediately with a plenty of water and seek medical attention.
 - Use only the recommended batteries. Do not mix old and new batteries. Always replace weak or empty batteries with a new, complete set of batteries at full capacity. Do not use batteries from different brands, types or with different capacities. Remove the batteries from the unit if it is not to be used for a long period of time, or if the batteries are empty!
 - When inserting the batteries, ensure that the polarity (+/-) is correct.
 - Never recharge normal, non-rechargeable batteries. This could lead to explosion during the charging process.
 - Rechargeable batteries are only to be charged under adult supervision.
 - Rechargeable batteries are to be removed from the toy before being charged.
 - The terminals must not be short-circuited.
 - Do not disassemble the device. In the event of a defect, please contact your dealer. The dealer will contact the Service Centre and can send the device in to be repaired, if necessary.
 - Tools with sharp edges are often used when working with this device. Because there is a

risk of injury from such tools, store this device and all tools and accessories in a location that is out of the reach of children.

- Keep instructions and packaging as they contain important information.

Parts overview (Fig. 1+2):

- ❶ 5x WF Eyepiece
- ❷ 16x WF Eyepiece
- ❸ Barlow lens
- ❹ Eyepiece socket
- ❺ Microscope head
- ❻ Objective-revolver
- ❼ Microscope stage
- ❽ Focus wheel
- ❾ LED lighting (transmitted light)
- ❿ On/off switch (illumination)
- ⓫ Microscope base
- ⓬ Battery compartment (3x AA)
- ⓭ 5 Slides and 5 preparations
- ⓮ Dimmer
- ⓯ Condenser
- ⓰ Adjusting lever for iris diaphragm
- ⓱ Microscope tools
- ⓲ Petri dish
- ⓳ Empty bottles
- ⓴ Cover plates
- ⓵ Smartphone holder

1. General/Location

First you must make sure that your microscope is on a stable and solid surface.

2. Electric LED lighting with dimmer

In the base of the microscope there is a battery compartment (12). Loosen the screws at the battery compartment cover with a small screw-

driver and remove the cover.

Place the batteries (3x AA) in the compartment so that the flat minus poles (-) press against the spring terminal and the plus poles (+) are touching the flat contact sheets.

Close the battery compartment with the cover and turn the microscope around again.

As your device has continuously controllable lighting (dimmer) optimal illumination of the object to be viewed is guaranteed.

3. Microscope setup

The microscope's eyepiece will now be prepared for the first observation.

First, rotate the eyepiece into a convenient position. Begin every observation with the lowest magnification. Place the microscope stage (Fig. 1, 7) with the focus wheel (Fig 1, 8) into the lowest position and rotate the objective revolver (Fig. 1, 6) until it locks on the lowest magnification (4x).

NOTE:

First, drive the Microscope's table (Fig. 1, 7) in its lowest position before changing the objective in order to prevent damages.

Insert the 5x eyepiece (Fig. 3, 1) in the Barlow lens (Fig. 3, 3).

Take care, that the Barlow lens is inserted completely in the eyepiece socket (Fig. 3, 4).

4. Condenser

The condenser is used to adjust the light bundling to the observed object in an optimal way. The iris diaphragm is used to increase the depth of focus. With an open aperture more light is getting through but the depth of focus is less. With partially closed aperture the depth of focus is higher but it is necessary to increase the brightness of the illumination then. It is recommended not to close the iris diaphragm too much when working with higher magnifications. Otherwise the image resolution will decrease. Move the lever for the iris diaphragm to right or left to change the aperture setting.

5. Observation

After you have set up the microscope with the corresponding illumination, the following principles are important:

Begin each observation with a simple observation at lowest magnification, so that the centre and position of the object to be viewed is in focus.

The higher the magnification the more light is required for good picture quality.

Place a permanent slide culture (Fig. 4) directly under the microscope lens on the plate (Fig. 4, 7) and clamp it tight. The specimen to be examined must be directly over the lighting.

Look through the eyepiece (Fig. 1, 1/2) and turn carefully the focusing wheel (Fig. 1, 8) until you can see a sharp picture.

Now you can get a higher magnification, while you pull out slowly the Barlow lens (Fig. 5, 3) of the eyepiece socket (Fig. 5, 4). With nearly entirely pulled out Barlow lens the magnification is raised to 2x.

For still higher magnification you can put the 16x eyepiece (Fig. 5, 2) into the objective revolver (Fig. 1, 6) and set on higher position (10x / 40x).

TIP:

Depending on the preparation higher magnifications do not always lead to better pictures.

With changing magnification (eyepiece or objective lens changes, pulling out of the Barlow lens) the sharpness of the image must be newly defined by turning the focusing wheel (Fig. 1, 8).

NOTE:

Please be very careful when doing this. When you move the mechanical plate upwards to fast the objective lens and the slide can touch and become damaged.

6. Viewed Object – condition and preparation

6.1. Condition

With the Barlow lens nearly fully extended magnification can be doubled. Both transparent and non-transparent specimens can be examined with this microscope, which is a transmitted light model. If translucent (transparent) specimens are examined the light from below goes through the specimen, lens and eyepiece to the eye and is magnified en route. Many small organisms of the water, plant parts and finest animal components have now from nature these transparent characteristic, other ones must be accordingly prepared. Is it that we make it by means of a pre-treatment or penetration with suitable materials (media) transparent or thus that we cut finest wafers off of them (hand cut, MicroCut) and these then examine. With these methods will us the following part make familiar.

6.2. Creation of thin preparation cuts

Specimens should be sliced as thin as possible, as stated before. A little wax or paraffin is needed to achieve the best results. A candle can be used for the purpose. The wax is put in a bowl and heated over a flame. The specimen is then dipped several times in the liquid wax. The wax is finally allowed to harden. Use a MicroCut or knife/scalpel (carefully) to make very thin slices of the object in its wax casing.

DANGER!

Be extremely careful when using the knives/scalpels or the MicroCut. There is an increased risk of injury due to the sharp edges!

These slices are then laid on a glass slide and covered with another.

6.3. Creation of an own preparation

Put the object which shall be observed on a glass slide and give with a pipette (Fig. 6, 17a) a drop of distilled water on the object (Fig. 6).

Set a cover glass (in each well sorted hobby shop available) perpendicularly at the edge of the water drop, so that the water runs along the cover glass edge (Fig. 7). Lower now the cover glass slowly over the water drop.

7. Experiments

Use the following web link to find interesting experiments you can try out. How do you hatch brine shrimp? How to make salt crystals and much more. <http://www.bresser.de/downloads>

Microscope accessories can be found here: <http://www.bresser.de/P8859480>

Smartphone holder

1. Place the smartphone with the front camera lens centred over the adapter opening.
2. Press the smartphone case firmly onto the suction cups of the mounting plate.
3. Tighten the safety strap around the smartphone and the camera plate.
4. Switch on the smartphone and start the camera app. The image can be seen on the display through the lens.

5. If necessary, adjust the sharpness of the image using the focus setting of the lens and the rotary wheel for adjusting the height of the mounting plate.
6. Use the zoom function of the smartphone to show the object of observation in full format on the display.

NOTE:

Make sure that the smartphone cannot slip off the cradle. Bresser GmbH accepts no liability for damage caused by a dropped smartphone!

Care and maintenance

Separate the device from the power supply before cleaning (remove batteries)!

Only use a dry cloth to clean the exterior of the device.

NOTE:
Do not use any cleaning fluid to avoid damaging the electronics.

Clean the lenses (eyepiece and/or lens) with a soft and lint-free cloth only (e.g. microfibre).

NOTE:
Do not apply excess pressure to the cloth so as to avoid scratching the lenses.

To remove more stubborn dirt, moisten the cleaning cloth with an eyeglass-cleaning solution and wipe the lenses gently. Protect the device from dust and moisture! After use - in particular in situations of high humidity - let the device acclimatize for a short period of time, so that the residual moisture can dissipate.

Troubleshooting

Error

No picture visible

Solution

- switch light on
- put condenser lens in place
- adjust focus

Specifications

Magnification table

<i>Eye p.</i>	<i>Obj ect.</i>	<i>Magn.</i>	<i>with Barlow I.</i>
5x	4x	20x	40x
5x	10x	50x	100x
5x	40x	200x	400x
16x	4x	64x	128x
16x	10x	160x	320x
16x	40x	640x	1280x

EC Declaration of Conformity

CE Bresser GmbH has issued a „Declaration of Conformity“ in accordance with applicable guidelines and corresponding standards. The full text of the EU declaration of conformity is available at the following internet address:
www.bresser.de/download/9810302/CE/9810302_CE.pdf

Disposal

Dispose of the packaging materials properly, according to their type (paper, cardboard, etc). Contact your local waste disposal service or environmental authority for information on the proper disposal.

Do not dispose of electronic devices in the household garbage!


As per the Directive 2012/19/EC of the European Parliament on waste electrical and electronic equipment and its adaptation into German law, used electronic devices must be collected separately and recycled in an environmentally friendly manner.

In accordance with the regulations concerning batteries and rechargeable batteries, disposing of them in the normal household waste is explicitly forbidden. Please make

sure to dispose of your used batteries as required by law – at a local collection point or in the retail market. Disposal in domestic waste violates the Battery Directive.

Batteries that contain toxins are marked with a sign and a chemical symbol.


- ¹ battery contains cadmium
- ² battery contains mercury
- ³ battery contains lead

Warranty & Service

The regular warranty period is 2 years and begins on the day of purchase. To benefit from an extended voluntary warranty period as stated on the gift box, registration on our website is required.

You can consult the full warranty terms as well as information on extending the warranty period and details of our services at www.bresser.de/warranty_terms.

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© I. Schmitt-Menzel
WDR mediagroup GmbH
Bresser UK Ltd.