MICROCAM MINI FULL HD HDMI

Art. No. 5914170

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
We wish you every success with your new BRESSER MicroCam Mini Full HD HDMI.

Special features of the BRESSER MicroCam Mini Full HD HDMI:

- Sony IMX307(C) sensor for crystal-clear images
- 1920x1080 pixels (Full HD), 1/2.8" sensor (5.57x3.13)
- Image output via HDMI
- Integrated on-screen menu for controlling the camera
- Powerful ISP (Image Signal Processor) and other related processing functions
- C-mount thread for easy adaptability
- Easy setup - in addition to the microscope, all you require is an HDMI monitor, the supplied USB mouse and the on-screen menu in the camera.
- **Scope of delivery:** MicroCam, AC adapter (Output: DC 12V,1A), USB mouse, HDMI-compatible cable.

We recommend that you read this manual carefully and keep it handy while you work with the camera. It contains all relevant information on operation, user safety and maintenance. By following these instructions, you will ensure that your camera functions reliably and trouble-free for years to come.

Applications

The MicroCam MINI Full HD HDMI is designed to capture digital images from stereo and biological microscopes with an HDMI interface. It is suitable for the following applications:

- Scientific research, education (teaching, demonstration and academic exchange)
- Research, Digital Laboratory
- Industrial visual inspection (PCB inspection, IC quality control)
- Histology
- Food safety studies (observation and counting of microbial colonies)
- Aerospace

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For information on liability claims or service requests, see "Warranty Terms" and "Contacting Customer Service" in this manual. Errors and omissions excepted - technical data is subject to change.

Copyright: 2023 Bresser GmbH
 All rights reserved.

The reproduction of this documentation - including in part - in any form (e.g. photocopy, print, etc.) as well as the use

and distribution by means of electronic systems (e.g. image file, website, etc.) without the prior written permission of the manufacturer or the distributor, Bresser GmbH, is not permitted. It is also prohibited to use or transmit this documentation with electronic systems (e.g. image file, website, etc.). The designations and brand names of the respective companies used in this documentation are generally protected by trade, trademark and/or patent law in Germany, the European Union and/or other countries.

Validity note

This documentation is valid for the products with the following listed
 Article number: 5914170
 Instructions version: v082023a
 Manual designation: Manual_5914170_Micro-Cam_de-en_BRESSER_v082023a
 For inquiries, please provide this information.

1.2 About these instructions/ Intended use

About this instruction manual

These operating instructions should be regarded as part of the camera!

Read the operating instructions carefully before using the device. Keep this instruction manual in a safe place for future reference. If the device is sold or passed on, the operating instructions must be handed to the subsequent owner/user of the product.

Intended use

- This camera is intended for private and business use.
- It was developed for capturing digital images from microscopes.
- The camera is intended for indoor use only.

WARNING:

This device is not intended to be used by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and/or knowledge, unless they have been given supervision or instruction on how to use the device by a person responsible for their safety.

1.3 Checking the package contents

- Open the packaging with care.
- Avoid getting fingerprints and dust on the camera's blocking filter. These can severely affect the image quality.
- Avoid dropping and damaging the camera and accessories.
- Remove all components from the packaging and check for completeness according to the delivery note.
- Handle the camera and its accessories with care and protect them from strong impacts and vibrations at all times.
- Keep the original packaging for storage purposes or returning the camera to the supplier in case repair or inspection is required.

1.4 Safety Instructions

DANGER OF AN ELECTRIC SHOCK!

This camera is powered by an AC adapter. If used improperly, there is a risk of electric shock. Electric shock can cause serious or fatal injuries. It is therefore essential to observe the following safety information:

- Disconnect the power supply when you are not using the camera, during long breaks, and before all maintenance and cleaning work (disconnect the mains plug).
- Place your camera in such a way that it can be disconnected from the power supply at any time. The mains socket should always be close to the camera and easily accessible to enable quick disconnection from the mains.

- To disconnect the camera from the mains supply, always pull on the mains plug and never on the cable!
- Check the camera, cables and connections for damage before use. Never use a damaged camera or power adapter.
- Damaged parts must be replaced immediately by an authorized service centre.
- Only operate the camera in a completely dry environment.
- The camera is equipped with a power pack (output: 12V,1A)
- To avoid electric shock, connect the supplied power cord to a properly grounded outlet.

DANGER OF SUFFOCATION!

Improper use of this product may result in suffocation, especially for children. It is therefore imperative that you observe the following safety information.

- Keep packaging materials (plastic bags, rubber bands, etc.) away from children! There is a danger of choking!
- Keep cables, straps and similar materials away from children - risk of suffocation!

EXPLOSION HAZARD!

Improper use of the camera and accessories can lead to an explosion!

It is essential that you observe the following safety information in order to avoid fires.

- Do not expose the camera and accessories to high temperatures and extreme temperature fluctuations.
- Only use the power pack supplied by BRESSER. Do not short-circuit the camera or throw it into a fire!
- Excessive heat and improper handling can cause a short circuit, fire or even an explosion!
- Do not use the camera and the supplied accessories in explosive environments, or in the presence of flammable solvents

such as alcohol, gasoline, or volatile anaesthetic agents.

NOTE: **DANGER OF MATERIAL DAMAGE!**

Improper handling may damage the camera and/or accessories. For this reason, only use the camera in accordance with the following safety information:

- Observe all safety instructions and the operating instructions for the devices to which the camera is connected.
- Use in commercial facilities: always comply with the applicable occupational health and safety regulations for electrical systems and equipment.
- Use in schools, training institutions, hobby or DIY workshops: Monitoring and operation exclusively by trained personnel.
- Do not disassemble the camera! Do not attempt to disassemble or replace connectors, optical systems or mechanical parts yourself. In the event of a defect, please contact your dealer. The dealer will contact the service centre and send the device for repair if necessary. The camera and its accessories should only be serviced and repaired by trained personnel.
- Do not expose the camera and its accessories to high temperatures, extreme temperature changes, high humidity, water or chemicals (acids and vapours, alkaline or caustic and corrosive materials).
- Protect the camera from grease and oils.
- Protect the camera from direct sunlight and direct proximity to radiators.
- Do not subject the camera to physical shock!
- Only use accessories and spare parts for the camera that correspond to the technical specifications.
- Always use the mains cable supplied by BRESSER. If an unsuitable mains cable is used, BRESSER cannot guarantee the

electrical functionality and safety of the camera.

- Only use the camera and the original accessories for the applications described in this manual.
- The manufacturer assumes no liability for any other application. This also applies to all repair and service work that is not carried out by authorised service personnel. Therefore all guarantee / warranty claims expire.

1.5 Cleaning and Maintenance

- The use of unsuitable cleaning agents and methods can damage the camera housing!
- Clean the outside of the camera housing with a dry cloth.
- Blow loose dust off the sensor surface.
- Use high quality lens cleaning towels or soft cloth and moisten it with a little bit of pure alcohol (available in a drug store). Gently wipe the sensor surface in a circular motion from the centre to the edge of the lens, applying a small amount of pressure.

1.6 Care, storage and transport

- Protect the camera from dust, dirt and moisture. This affects the quality of the results.
- Avoid fingerprints and similar contamination on the optical surface.
- Switch off the device after each session and cover it (to protect against dust and moisture).
- If you do not plan to use the camera for an extended period, store the camera and accessories in a closed container in a dry and mould-free location.
- If the camera and accessories are used in a warm and humid climate, special care is required to prevent the formation of mould.

- Use the original packaging for shipping or transport.

2. Scope of delivery

The package contents are provided for reference only and may vary.

Fig. 1: Package contents of the MicroCam MINI Full HD HDMI.

A: MicroCam MINI Full HD HDMI

B: HDMI compatible cable

C: USB mouse

D: AC adapter: Entrance: AC 100~240V 50Hz/60Hz, Output: DC 12V,1A

3. Product information

- Sony Starvis "back-illuminated" CMOS sensor (Back-illuminated: the sensor is illuminated from behind - the side opposite the conductor. The light rays hit the pixel unhindered, which increases the sensitivity to light)
- FHD HDMI-compatible video outputs
- SD card slot for storing the captured images and videos
- Integrated on-screen menu for controlling the camera
- Powerful ISP (Image Signal Processor) and other related processing functions.

Order code	MicroCam MINI Full HD HDMI
Sensor/ sensor size	Sensor: Sony IMX307(C) Sensor size (mm): 1/2.8" (5.57x3.13)
Pixel (µm)	2.9x2.9
G sensitivity Dark signal	1300mV with 1/30s
FPS/Resolution	60@1920*1080 (HDMI compatible)
Binning	1x1
Exposure time (ms)	0.01~1000

List of abbreviations:

Pixel: The resolution of a digital camera is given in pixels. A pixel is a single point in the image and contains information about colour and brightness

FPS: Frames per second, frame rate

Binning: The grouping together of adjacent image points (pixels) in a digital camera

Connections and functional description

Interface/ Display	Description
HDMI	Compliant with HDMI 1.4 standard, 1080P video output for standard FHD monitor
LED	LED status indicator, illuminates when the camera is in operation
USB mouse	Connection socket for USB mouse to control the on-screen menu
SD	Slot for an SD card for storing images and videos; conforms to the SDIO3.0 standard
DV12V	DC connection socket (12V,1A)


Fig. 2: Available connections on the back of the MicroCam MINI Full HD HDMI

Video output interface	Description
HDMI compatible interface	Complies with HDMI-compatible 1.4 standard; 60fps@1080P
Other functions	Description
Save video	Video format: 2M (1920*1080) H264-encoded MP4 file Video storage frame rate: 50-60fps (depends on SD card performance)
Image capture	2M (1920*1080) JPEG image on SD card
Save measurements	Measurement information is stored in different layers with image content. The measurement data is saved together with the image content in <i>burn-in mode</i> (<i>objects containing the measurements are integrated into the image, and cannot be edited by the user</i>).
ISP function (Image Signal Processor)	Exposure (automatic/manual exposure) Gain White balance (manual/auto/ROI mode) Sharpen 3D denoise (noise reduction) Saturation adjustment Contrast adjustment Brightness adjustment Gamma value adjustment (adjusts the contrast range of the midtones in the image) Colour image to grayscale image (Color to Gray) 50Hz/60Hz anti-flicker function
Image editing	Zoom in/out Mirror/flip Freeze Cross line Overlay Browser for saved files Video playback Measurement function
Restore factory settings	Resets the camera parameters to the factory settings
Support multiple languages	English/Simplified Chinese/Traditional Chinese/Korean/Thai/French/German/Japanese/Italian/Russian

Operating conditions	
Operating temperature	-10°C ~ 50°C
Storage temperature	-20°C ~ 60°C
Operating humidity	30~80%RH (Relative Humidity)
Storage humidity	10~60%RH (Relative Humidity)

4. Commissioning the camera

Easy setup - in addition to the microscope, all you require is an HDMI monitor, the supplied USB mouse and the built-in camera software.


Fig. 3: Connecting the MikroCam MINI Full HD HDMI to an HDMI monitor.

4.1 Step-by-step setup instructions:

A Attach the camera to the microscope

The MicroCam MINI Full HD HDMI features a C-Mount thread for connecting the microscope. This makes it compatible with various microscopes and adapters.

Bresser offers 2 different adapters:

- BRESSER Science C-Mount MicroCam Adapter (**Item No. 5942101**) - for connecting a BRESSER MicroCam or other C-Mount cameras to a wide range of BRESSER Science microscopes.
- BRESSER Photo Adapter 23.2 mm / C-Mount (**Item No. 5942070**) - This adapter is required for connecting a MicroCam or C-Mount camera if the microscope has a 23.2 mm tube.

The camera (1) is screwed directly into the C-mount adapter (2) on the microscope head.

NOTE:

The MicroCam MINI Full HD HDMI typically has a smaller field of view than the eyepieces.

Sensors smaller than 1" show a smaller field of view than a 10x eyepiece with a field number of 22.

Remember:

The microscope and the associated settings, including the KÖHLER illumination, should be configured before the camera is installed.

B Connect the camera to the monitor

Connect the camera (1) to an HDMI monitor (4) using the HDMI cable (3).

C Connect the camera to the USB mouse

Connect the supplied USB mouse (5) to the USB port on the camera (1).

Remember:

If the mouse is not connected, the camera is only operated in automatic mode. In this case, it is not possible to record images.

D Insert the SD card

Insert the supplied SD card (6) into the SD card slot on the camera (1).

Remember:

If no SD card is inserted, only the live image will be displayed, and it will not be possible to take pictures and videos.

E Power supply

Connect the supplied power adapter to the camera (1) using the DC plug (7), and then connect the power supply (8) to a mains socket. The LED status indicator (9) will turn on.

Turn on the monitor and look at the image on the screen.

5. User interface of the camera

Move the mouse left, up or down on the interface. The camera displays a range of menu options, which you can select with the mouse.

The camera user interface includes:

- The **CAMERA CONTROL PANEL (1)** on the left side of the screen
- The **UPPER MENU BAR (2)** (for taking measurements)
- The **BOTTOM TOOLBAR (3)** (for controlling the camera)

NOTE:

In addition, a ruler appears at the top and left edge of the screen, which can be moved.



Fig. 4: User interface of the MikroCam Mini Full HD HDMI.


5.1 Camera control panel

- The **camera control panel** allows you to configure the camera settings so that you can obtain the best image or video quality depending on the application.
- **It appears automatically when the mouse pointer is moved to the left side of the screen.**
- Left-click the 📌 icon to permanently show or hide the camera control panel.
- Sliders which appear blue can be moved with the mouse. To move a slider, left click the corresponding slider. and drag it to the desired position.

Remember:

The control panel of the camera is not automatically displayed during the measurement process. After you have completed a measurement, the upper menu bar can be closed using the **X** button, and the camera control panel can be displayed again on the left side of the screen.

Camera control panel	Function	Description
 <p>Camera Control Panel</p> <p><input checked="" type="checkbox"/> Auto Exposure:</p> <p>Exposure Compensation: 56</p> <p>Exposure Time: 1ms</p> <p>Gain: 90</p> <p>Red: 116</p> <p>Green: 102</p> <p>Blue: 79</p> <p><input type="radio"/> Auto <input checked="" type="radio"/> Manual <input type="radio"/> ROI</p> <p>Sharpness: 28</p> <p>Denoise: 8</p> <p>Saturation: 31</p> <p>Gamma: 9</p> <p>Contrast: 29</p> <p>Brightness: 48</p> <p><input checked="" type="radio"/> DC <input type="radio"/> AC(50Hz) <input type="radio"/> AC(60Hz)</p> <p>Default</p>	Single image	Capture image and save to SD card (note: an SD card must be inserted).
	Record	Record video and save to SD card.
	Automatic exposure	<i>If automatic exposure is enabled, the exposure time and sensitivity will be set automatically.</i>
	Exposure compensation	<i>Available when automatic exposure is enabled. Move the slider left or right to adjust the exposure compensation according to the current image brightness and set the appropriate brightness value (0-255).</i>
	Exposure time	<i>Available when automatic exposure is disabled. Move the slider left or right to shorten or lengthen the exposure time (0.01ms-1000ms) and adjust the brightness of the image.</i>
	Gain	<i>Available when automatic exposure is disabled. Adjust the gain to decrease or increase the brightness of the image (0-90). Noise is reduced or increased accordingly.</i>
	Red	Move the slider left or right to decrease or increase the amount of red in the image (0-400).
	Green	Move the slider left or right to decrease or increase the amount of green in the image (0-400).
	Blue	Move the slider left or right to decrease or increase the amount of blue in the image (0-400).
	Automatic white balance (Automatic)	Adjust the white balance according to the image.
Manual white balance	In addition to automatic white balance, the sliders for the red and blue channels can be adjusted manually.	
ROI white balance (ROI)	The white balance can be changed within the ROI (Region of Interest).	

Camera control panel	Function	Description
	Sharpness	Adjust the sharpness (0-100): (Note: a high value does not automatically lead to a better result.)
	Noise reduction	Adjust the noise reduction setting (0-50): In low light, this value can be increased to reduce image noise.
	Saturation	Adjust the saturation: use this option to increase or decrease the colour of the image (0-100).
	Gamma	Adjust the brightness values: this setting allows you to change the image brightness, primarily in the medium brightness ranges (1-20).
	Contrast	Adjust the contrast (0-100): (a high value does not automatically lead to a better result.)
	Brightness	Adjust the image brightness (0-100)
	DC/AC (50HZ)/ AC (60HZ)	This function enables you to suppress image flickering, depending on the type of lighting or mains frequency. <ul style="list-style-type: none"> • DC: Direct current • 50 50 Hz alternating current • 60:00:00 60 Hz alternating current
	Standard	Clicking on this button resets all the functions to automatic. This function is very useful if manual adjustments were unsuccessful and the picture is unsatisfactory.

5.2 Top menu bar (for taking measurements)

- The top toolbar appears when you move the mouse pointer near the top of the image window.
- Left-click the  icon to permanently show or hide the top menu bar.

Remember:

If the top menu bar is pinned to the screen, the camera control panel and bottom toolbar will be greyed out. To end the measurement process and access the other toolbars, click the  icon again or close the upper menu bar with the  icon.

- During the measurement process, you can select a specific measurement object. The object position and properties control bar  will be displayed, which allows you to change the object position and properties of the selected objects.

The following section shows the various functions of the top menu bar:

Icon	Name/function
	Pin (auto hide) Pin/hide the top menu bar (NOTE: other menu bars are greyed out)
	Visible Display/hide the measurement functions
	Selecting the unit of measurement for the measurements: <ul style="list-style-type: none"> • Pixel • Nanometres (nm) • Micrometres (µm) • Millimetres (mm) • Centimetres (cm) • Metres (m) • Inches (in) • Million (mils)
	Object selection Click the "arrow" to move or change measurements.
	Select magnification for measurement AFTER calibration: This value should always match the selected magnification, otherwise the measurement will be incorrect.
	Angle (measurement function)
	Four point angle (measurement function)
	Point (measurement function)

	Any distance (measurement function)
	Three-point distance (measurement function)
	Horizontal distance (measurement function)
	Vertical distance (measurement function)
	Parallel (measurement function)
	Three point vertical line (measurement function)
	Rectangle (measurement function)
	Rectangle (Three points) (measurement function)
	Ellipse (measurement function)
Icon	Name/function
	Five-point ellipse (measurement function)
	Centre circle+radius (circle starting from the centre) (measurement function)
	Circle (three points) (measurement function)
	Ring (of two circles) (measurement function)
	Ring (consists of two circles and three points) (measurement function)
	Two circles (and their centre distance) (measurement function)

Icon	Name/function
	Three points, two circles (Two circles, each with three points and their centre distance) (measurement function)
	Circular arc (measurement function)
	Text (insert text field) (additional function)
	Polygon (additional function)
	Curve (drag freely with the pointer) (additional function)
	Scale bar (additional function)
	Arrow (additional function)
	Calibration Perform a calibration to determine the appropriate relationship between magnification and resolution, and on this basis, calculate the relationship between the measurement unit and the pixel size of the sensor. Calibration must be performed using a micrometer.
	Export to CSV Save all measured values in a CSV file (*.csv)
	Setting Configure the settings for displaying lines and text.
	Delete all measurement objects All measurements are cleared from the screen. A prompt will be displayed in which you can confirm your selection.
	End Exit/hide measurement mode

Control bar for object position and properties

Opens **AFTER THE MEASUREMENT IS COMPLETED**

Left-click on a measurement object to display the control bar for the object position and properties. The user can move the object by dragging it with the mouse. More precise movements can be obtained using the control bar. To do this, press and hold the desired button.
Click on the 3 colour points to adjust the colour.

 to the left
 to the right
 upwards
 downwards
 colour matching
 delete

5.3 Bottom toolbar

- The **bottom toolbar** appears when you move the mouse pointer near the **bottom of the image window**.
- The toolbar is used to control the camera.

Icon	Function
	Zoom in
	Flip horizontally
	Colour conversion (colour/grey)
	Grid Display a crosshair
	Settings
	Zoom Out
	Flip vertically
	Freeze (Video freeze) Freezes the image

Icon	Function
	Browser Browse the images and videos on the SD card
	Info Version info

5.4 ✂ Settings button

The sub-menus in **Settings** are more complex than the other functions. The following section provides an overview of each sub-menu:

5.4.1 Settings > Measurement


Fig. 5. Structure of the measurement.

Global	Accuracy	Save calculation results to specified number of decimal places (0-5)
	Font size	Use this option to set the font size (small, medium, large)
	Cursor	Use this option to define the type of cursor (none/single cross) and the colour
	Miscellaneous	Hide labels when moving the measurement object.
Calibration	Line width	Used to specify the line width (1-10)
	Colour	Used to set the line colour. Click the coloured dots button to change the colour.
Angle, Point, Line, Parallel, Vertical Line, Rectangle, Ellipse, Circle, Ring, Two Circles, Arc, Text, Polygon, Curve, Scalebar, Arrow		Select one of the listed commands with the mouse and click on the + icon. The corresponding properties of the command are thus displayed. (Label type: Used to define the shape more precisely)

5.4.2 Settings > Magnification


Fig. 6: Magnification settings.

Name	The magnification name, typically the magnification of the microscope objective, is used as the magnification name during calibration (e.g. 4x, 10x, 100x, etc). Other custom information can also be added to the magnification name, e.g. the microscope model, the name of the user, etc.
Disbandment	Pixels per millimetre (mm). Imaging devices such as microscopes have a high resolution value.
Clear all	Clear all calibrated magnifications.
Delete	Delete the selected magnification.

5.4.3 Settings > Image format


Fig. 7: Image format settings.

Storage method for measurement objects	<p>Fusion mode: Measuring objects are integrated into the image, and cannot be edited by the user.</p> <p>Layer mode: Measurement objects are saved in different layers with image data in the target file. The user can edit the measurement objects in the target file with the software on the PC. This mode is reversible.</p>
---	--

5.4.4 Settings > Video


Fig. 8a: Video settings: Video playback.

Video playback	Fast forward/rewind The time interval (0 to 600 seconds) for playing video files.
-----------------------	---


Fig. 8b: Video settings: Video encoding.

<p>Video encoding</p>	<p>H264: The encoding format of the video files is H264 format.</p> <p>H265: The encoding format of the video files is the H265 format (requires hardware support). Check beforehand whether your playback device supports this standard.</p>
------------------------------	---

5.4.5 Settings > Storage

Fig. 9 File system format of the storage device.

<p>Storage device</p>	<p>SD card: Only the SD card is supported as a storage device</p>
<p>File system format of the storage device</p>	<p>Lists the file system format of the current storage medium: FAT32: The SD card file system is FAT32. The maximum video file size of a single file is 4G bytes exFAT: The SD card file system is exFAT. The maximum video file size of a single file is 4G bytes NTFS: The SD card file system is NTFS. The maximum video file size of a single file is 4G bytes. Use a PC to format the SD cards and switch between FAT32, exFAT and NTFS. Unknown status: SD card is not recognized or the file system is not identified.</p>

5.4.6 Settings > Files


Figure 10: Image file settings.

Image file name	<p>Automatically: The image files are automatically saved with the specified prefix.</p> <p>Manual: The user must specify the filename BEFORE saving the image.</p>
Video file name	<p>Automatically: The video file is automatically saved with the specified prefix.</p> <p>Manual: The user must specify the name of the video file BEFORE recording the video.</p>
<p>NOTE: The maximum video file size is 4G bytes. In the case of extended video recordings, several video files can be generated automatically.</p>	

5.4.7 Settings > Language

Fig. 11 Language selection.

English (English)	Set the language of all software to English.
Simplified Chinese	Set the software language to Simplified Chinese.
Traditional Chinese (Traditional Chinese)	Set the software language to Traditional Chinese.
Korean (Korean)	Set the software language to Korean.
Thai (Thai)	Set the language of all software to Thai.
French (French)	Set the software language to French.
German (German)	Set the software language to German.
Japanese (Japanese)	Set the software language to Japanese.
Italian (Italian)	Set the software language to Italian.
Russian (Russian)	Set the software language to Russian.

5.4.8 Settings > Miscellaneous


Fig. 12 Other setting options.

Ruler	Choose whether to display the ruler on the screen. When the measurement is calibrated, the ruler changes according to the number of pixels.
Automatic exposure	The maximum exposure time during auto exposure can be set (1-100). Setting this option to a lower value can ensure a faster frame rate during auto exposure.
ROI	Select the colour of the ROI - Click on the button with the three coloured circles
Camera parameters	Import: Import camera parameters from the SD card to use the previously exported camera parameters. Export: Export camera parameters to the SD card to use the previously exported camera parameters.
Reset to factory settings	Resets the camera parameters to the factory settings

6. Product dimensions

The camera

The housing of the MicroCam MINI Full HD HDMI is made from a durable CNC-milled aluminium alloy. A high-quality IR-CUT filter is used to protect the camera sensor. The IR-CUT filter only lets visible light through and unwanted IR light is reflected. The camera contains NO moving parts, making it robust and durable.


Fig. 13:00 Product dimensions (mm) of the MikroCam MINI Full HD HDMI

1" (32): The C-mount connection thread has an actual diameter of 1 inch (25.4 mm) and a thread pitch of 1/32 inch.

¼" (20): Tripod connection thread (camera tripod).

7. EC declaration of conformity

 Bresser GmbH has prepared a “CE Declaration of Conformity” in accordance with the applicable guidelines and corresponding standards. The full text of the CE declaration of conformity is available at the following Internet address:
www.bresser.de/download/5201000/CE/5201000_CE.pdf

8. Disposal

Dispose of the packaging materials properly, according to their type, such as paper or cardboard. Contact your local recycling service or environmental authority for information about the proper disposal.

Do not dispose of electronic devices in the household garbage! According to the European Directive 2002/96/EC on Waste Electrical and Electronic Equipment and its transposition into national law, used electrical equipment must be collected separately and recycled in an environmentally sound manner.

9. Warranty & Service

The regular guarantee period is 2 years and begins on the day of purchase. To benefit from an extended voluntary guarantee period as stated on the gift box, registration on our website is required.

You can consult the full guarantee terms as well as information on extending the guarantee period and details of our services at www.bresser.de/garantiebedingungen.

10. UKCA DECLARATION OF CONFORMITY

A “Declaration of conformity” in accordance with the applicable directives and corresponding standards has been prepared by Bresser GmbH. This can be viewed at any time upon request at www.bresser.de/download/5914170/UKCA/5914170_UKCA.pdf

BRESSER UK Ltd. • Suite 3G, Eden House
Enterprise Way, Edenbridge, Kent TN8 6HF, GB

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GB IE

Please contact the service centre first for any questions regarding the product or claims, preferably by e-mail.

E-Mail: service@bresseruk.com
Telephone*: +44 1342 837 098

BRESSER UK Ltd.

Suite 3G, Eden House
Enterprise Way
Edenbridge, Kent TN8 6HF
Great Britain

*Number charged at local rates in the UK (the amount you will be charged per phone call will depend on the tariff of your phone provider); calls from abroad will involve higher costs.

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Bresser UK Ltd.